

Revision 141 software requires these hardware modifications for the MOTOR BREAK IN and MOTOR CURRENT test to function properly:

